



# MEDSTONE

# Zogear<sup>®</sup> Surgical Face Mask Type I



Mask is fluid resistant and provides the wearer protection against large droplets, splashes, or sprays of bodily or other hazardous fluids. Protects others from the wearer's respiratory emissions.

## Specifications:

- BFE:  $\geq 95\%$ \*
- Breathability (Delta P):  $< 4.0 \text{ H}_2\text{O}/\text{cm}^2$
- Flammability: Class I

## Characteristics:

- 100% Latex-free soft stretchable ear-loops
- Fluid resistant white inner layer
- Adjustable nose piece
- High-standard breathability, fiberglass free, hypoallergenic
- Comfortable inner lining

## Certifications:

- CE
- FDA
- ISO13485

\*Type I based on EN14683



MEDSTONE

Medstone  
Amstelplein 28A  
Mondriaan Tower (13th Floor)  
1096 BC Amsterdam  
The Netherlands

Company website: [www.medstone.nl](http://www.medstone.nl)  
Online shop: [shop.medstone.nl](http://shop.medstone.nl)  
Contact: [info@medstone.nl](mailto:info@medstone.nl)

Disclaimer: No rights can be derived from the information provided on this sheet  
©2020 All Rights Reserved.  
Medstone<sup>®</sup> is a registered trade name of Fresco Farma B.V.