



# MEDSTONE

# Zogear<sup>®</sup> Surgical Face Mask



Mask is fluid resistant and provides the wearer protection against large droplets, splashes, or sprays of bodily or other hazardous fluids. Protects others from the wearer's respiratory emissions.

## Specifications:

- Fluid/splash resistance pressure: 160mmHg/21.3 kPa\*
- PFE: ≥ 98%
- BFE: ≥ 98%\*
- Breathability (Delta P): < 5.0 H<sub>2</sub>O/cm<sup>2</sup>
- Flammability: Class I

## Characteristics:

- 100% Latex-free soft stretchable ear-loops
- Fluid resistant white inner layer
- Adjustable nose piece
- High-standard breathability, fiberglass free, hypoallergenic
- Comfortable inner lining

## Certifications:

- CE
- FDA
- ISO13485

\*Type IIR based on EN14683 / Level 3 based on ASTM F2100 testing standards



MEDSTONE

Medstone  
Amstelplein 28A  
Mondriaan Tower (13th Floor)  
1096 BC Amsterdam  
The Netherlands

Company website: [www.medstone.nl](http://www.medstone.nl)  
Online shop: [shop.medstone.nl](http://shop.medstone.nl)  
Contact: [info@medstone.nl](mailto:info@medstone.nl)

Disclaimer: No rights can be derived from the information provided on this sheet  
©2021 All Rights Reserved.  
Medstone<sup>®</sup> is a registered trade name of Fresco Farma B.V.